

BL224 Touch

Electronic Analytical Balance

USER MANUAL

\boldsymbol{I} About the Balance

Specification

Model	Range	Division	Repeatability	Linearity	Pan Size
BL224	220g	0.1mg	0.2mg	0.5mg	ø 80mm
Touch					

Features

Stainless Steel Pan

Super 5 inch Touch panel display

RS232 Interface

Mains adapter supplied as standard

Windshield supplied as standard

Height adjustable feet

Internal Auto Calibration

Selectable measure units:mg, g, oz, ct


Memory for accumulated time

Applications

Weighing Net weight / tare Below balance weighing Piece counting function Density Test Percentage Test

9.1 DATA OUTPUT

1	Model or a decimal point				
2	A space or a decimal point				
3	A space or *				
4	+ or - or a decimal point				
5	data				
6	Data or a decimal point				
7	Data or a decimal point				
8	Data or a decimal point				
9	Data or a decimal point				
10	Data or a decimal point				
11	Data or a decimal point				
12	Data				
13	Unit 1				
14	Unit 2				
15	Unit 3				
16	Enter				
17	Wrap				



We can choose the suitable baud rate and communication style for different print or connect to the computer

8.1 FAULT JUDGMENT

reasons.

If you come across any problem, you can check the

RESON	EXCLUDE
Not connected to the power	Plug in the power line;
supply;	Replace the fuse;
Fuse is broken;	Replacement of power transformer;
Power supply damaged;	Contact factory
Bad working conditions;	Keep the environment stable, close
The wind screen is open;	the windows and doors;
Something between the	Close the glass door;
table and balance;	
The power unstable;	
Weighing unstable;	
The balance not calibrated.	Calibration.
No tare before weighing.	Tare before weighing.
No well leveled.	Adjust the level feet.
	supply; Fuse is broken; Power supply damaged; Bad working conditions; The wind screen is open; Something between the table and balance; The power unstable; Weighing unstable; Weighing unstable; No tare before weighing.

► 10

weigihing

IISETTINGS

TABLE OF CONTENTS

1.0 KNOW YOUR BALANCE 2.0 SPECIFICATIONS **3.0 UNPACKING THE BALANCE 4.0 CALIBRATION 5.0 LOCATING THE BALANCE 6.0 SETTING UP THE BALANCE 6.1 ASSEMBLING THE BALANCE 6.2 LEVELLING THE BALANCE** 7.0 OPERATION SETTINGS 7.1 SENSITIVITY SETTING **7.2 SPEED SETTING 7.3 WEIGHING UNIT SETTING 7.4 COUNTING FUNCTION** 7.5 PERCENTAGE FUNCTION 7.6 DENSITY FUNCTION **8.1 FAULT JUDGMENT** 9.1 DATA OUTPUT



1.0 KNOW YOUR BALANCE

Thank you for selecting the BL 224 Touch Analytical Balance. This Instruction Manual will guide you of the installation, accessories, trouble-shooting, after sales service information, general maintenance of the balance, etc. it will also guide you through the various applications. Please read this Manual thoroughly before starting the operations. If you need any clarifications, feel free to contact us.


The BL 224 Touch balances are ideal for laboratory and general purpose weighing. The balances can also be used for some advanced weighing functions.



FEATURES:

Stainless Steel Pan Super bright LCD display with backlight RS232 Interface Adapter supplied as standard Windshield supplied as standard Height adjustable feet, Bubble Level Internal Auto Calibration Selectable measure units, g, oz, ct,mg Memory for accumulated time Weighing function Piece counting function Below balance weighing

7. 5 PECENTAGE FUNCTION



1. Choose the 100% which you will put on the pan


2.Touch the weighing value you want weighing, here is the sample weight, after the weighing stable touch again, the result will check and store.

3.Touch here and weighing the total weight, and touch again check the weight.


4. The percentage will be show here.

7.6 Density Function.

Same operation with density function, just touch the air weight and weighing in the air, then weighing in the liquid, the result will show in the display.



01




7.3 WEIGHING UNIT SETTING

Under the weighing state, press "CONVERSION" to choose the unit. There are 4 units :g/ct/oz/mg.

7.4 COUNTING FUNCTION

There are counting functions in the balance, And please keep the products have the same weight, and the minimum weight must be ≥ 0.5 mg.



- 1. Choose the sample quantity which you will put on the pan
- 2. Touch the weighing value you want weighing, here is the sample weight, after the weighing stable touch again, the result will check and store.
- 3.Touch here and weighing the total weight, and touch again check the weight.
- 4. The quantity will be show here.

2.0 SPECIFICATIONS

Model	BL224 Touch		
Range	0-220g		
Division	0.1mg		
Repeatability	0.2mg		
Linearity	0.3mg		
Resp. Time	≤ 4 Sec.		
Pan Size	ø80mm		
Interface	RS-232 bi-directional		
Overall	340 x 215 x 350mm		
Dimensions			
Operating	+10°C to 40°C		
temperature	$+10^{\circ}$ C to 40° C		
Power Supply	Output voltage 5 VDC, 600mA through		
	External Power Adapter as standard (Input Voltage220 V)		
Net Weight	8.5 kg		

05

3.0 UNPACKING THE BALANCE

Remove the balance from the packing by carefully lifting it out of the box. Inside the box you will find everything needed to start using the balance

4.0 CALIBRATION THE BALANCE

Press"CAL" key, and it will show and shining: CAL-XXX.XXXX, and the balance begin internal calibration.

You also can choose manual calibration with standard weight, and there will have the notes tell you step by step:







One-point calibration Multi-Calibration In

Internal-Calibration

5.0 LOCATING THE BALANCE

• The balance should not be placed in a location that will reduce the accuracy.

• Avoid extremes of temperature. Do not place in direct sunlight or near air conditioning vents.

• Avoid unsuitable tables. The table or floor must be rigid ,no vibration.

• Avoid unstable power sources. Do not use near large users of electricity such as welding equipment or large motors.

- Do not place near vibrating machinery.
- Avoid high humidity that might cause condensation. Avoid direct contact with water. Do not spray or immerse the balances in water.
- Avoid air movement such as from fans or opening doors. Do not place near open windows or air-conditioning vents.

• Keep the balance clean. Do not stack material on the balances when they are not in use.

6.0 SETTING UP THE BALANCE


6.1 ASSEMBLING THE BALANCE

- Locate balance on solid surface, free from vibration
- Open the sliding door and gently place the stainless steel top.
- Level balance using the adjustable feet and the bubble level
- Connect the power to the balance
- For best performance, let the balance warm up for 30-60 min. and calibrate before using

6.2 LEVELLING THE BALANCE

After placing the balance in a suitable place, level it by using the bubble level. To level the balance turn the two adjustable feet at the rear of the balance until the bubble in the bubble level is center.

7.0 OPERATION SETTINGS



7.1 SENSITIVITY SETTING

