

pH AND CONDUCTIVITY BUFFER SOLUTION



makes pH and Conductivity easy

A wide range of solutions for pH, conductivity and maintenance is available in the form of practical 25ml sachets.

Single-dose sachets with long life time are quick and easy to use. Each sealed, light-tight sachet holds just right amount of solution. Every time your instrument is calibrated using these sachets, it like using a newly-opened bottle of solution.

Solutions in sachets are also supply with NIST certificate traceability.

All calibration solutions are provided with a label presenting a reference table of the relationship between pH or conductivity value and temperature, nominal value, batch of production and expire date.

Sachets solution are also available in one value carton box, 20 sachets, and combined kits for practicality. These kits are useful for multiparameter instruments or for two or three points calibrations and maintenance.



p/n	Description
-----	-------------

Single buffer

51102013	1x20 SACHET 25ml pH 4.01±0.01/25°C (red), with Certificate of Analysis
51102023	1x20 SACHET 25ml pH 7.00±0.01/25°C (red), with Certificate of Analysis
51102043	1x20 SACHET 25ml pH 10.01±0.02/25°C, with Certificate of Analysis
51102053	1x20 SACHET 25ml EC 1413 µS/cm ±1%/25°C, with Certificate of Analysis
51102063	1x20 SACHET 25ml EC 12880 µS/cm ±1%/25°C, with Certificate of Analysis
51102083	1x20 SACHET 10ml STORAGE Solution for pH/ORP sensors
51102093	1x20 SACHET 25ml Distilled Water

Calibration Set

51102203	Set CAL pH BIS 1x20 SACHET 25ml (10 x pH 4.01/7.00)
51102213	Set CAL pH TRIS 1x20 SACHET 25ml (7 x pH 4.01/7.00 + 6 x pH 10.01)
51102223	Set CAL COND 1x20 SACHET 25ml (10 x EC 1413/12880µS/cm)
51102233	Set CAL PC 1x20 SACHET 25ml (7 x pH 4.01/7.00 + 6 x EC 1413µS/cm)
51102243	Set REFILL pack for pH 1x20 SACHET 25ml (6 x pH 4.01/7.01/H2O + 2 x Storage 10ml)
51102253	Set REFILL pack for COND 1x20 SACHET 25ml (7 x EC 1413/12880µS/cm + 6 x H2O)
51102263	Set REFILL pack for PC 1x20 SACHET 25ml (5 x pH 4.01/7.01/EC 1413µS/cm + 3 x H2O + 2 x Storage 10ml)

